



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

ASPRV1, 191-326 aa, Human, His tag, E.coli

产品货号: TP01254

第三版

别名:Aspartic peptidase, retroviral-like 1, MUNO, SASP, SASPase, Taps.

描述:ASPRV1, also known as MUNO, SASP and ASAPase, is a protein that contains 1 peptidase A2 domain. This protein is expressed primarily in the granular layer of the epidermis and inner root sheath of hair follicles and localized to membrane region.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% glycerol, 1mM DTT

分子量:17.2 kDa (159aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSSMGKGYLKGKIGKVPVRFLVDSGAQVSVVHPNLWEEVTDGDLDLTLQPFENVVK
VANGAEMKILGVWDTAVSLGKLLKAQFLVANASAEAAIIGTDVLQDHNAILDFEHRCTLKGKKFRLLPVGGSLDEFD
LE

纯度:> 95% by HPLC

浓度:0.25 mg/ml

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.