



本公司提供的电版本说明书仅供参考，实验请以收到的纸质手册为准。

ATOX1, 1-68aa, Human, His tag, E.coli

产品货号: TP01267

第三版

**别名:**Copper transport protein ATOX1, Antioxidant protein 1, ATX1, HAH1.

**描述:**Copper transport protein ATOX1, also known as HAH1, a metal transport protein that belongs to the ATX1 family. ATOX1 is involved in cellular antioxidant defense and can bind and deliver cytosolic copper to the copper ATPase proteins.

**配方:**Liquid. In 20mM Tris-HCl buffer (pH8.0) containing 1mM DTT, 10% glycerol

**分子量:**9.5 kDa (88aa)

**序列:**

MGSSHHHHHHSSGLVPRGSHMPKHEFSVDMTCGGCAEAVSRVLNKLGGVKYDIDLPNKKVCIESEHSMDTLLATLKKT  
GKTVSYLGLE

**纯度:**> 95% by HPLC

**浓度:**1 mg/ml (determined by Bradford assay)

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.