

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

BAALC, 1-145aa, Human, His tag, E.coli

产品货号: TP01288

第三版

别名:Brain and acute leukemia, cytoplasmic isoform 1.

描述:Brain and acute leukemia, cytoplasmic (BAALC) is localizes to both the membrane and the cytoplasm and exists as multiple alternatively spliced isoforms. Expressed by hematopoetic and neural cells, BAALC interacts with CaMKII and is thought to play a role in synaptic function at postsynaptic lipid rafts.

配方:Liquid. In 20mM Tris-HCl buffer (pH8.0) containing 1mM DTT, 20% glycerol, 2mM EDTA

分子量:17.7 kDa (165aa)

序列:

MGSSHHHHHHSSGLVPRGSHMGCGRADAIIEPRYYESWTRETESTWLTYTSDAPPSAAAPDSGPEAGGLHSGMLEDG
LPSNGVPRSTAPGGIPNPEKKTNCETQCPNPQSLSSGPLTQKQNGLQTTEAKRDAKRMPEKVTINVTDSIQQMDRSTRIT
KNCVN

纯度:> 95% by HPLC

浓度:0.5 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.