

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

BATF, 1-125aa, Human, His tag, E.coli

产品货号: TP01297

第三版

别名:Basic leucine zipper transcriptional factor ATF-like, B-ATF, BATF1, SFA-2, SFA2.

描述:BATF, also known as SFA2, is a nuclear basic leucine zipper protein that belongs to the AP-1/ATF superfamily of transcription factors. BATF is strongly expressed in mature T and B lymphocytes, and is upregulated after transformation by human T-cell leukemia virus type I.

配方:Liquid. In 20 mM Tris-HCl buffer (pH8.0) containing 0.2M NaCl, 2mM DTT, 40% glycerol

分子量:16.2 kDa (145aa), confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMPHSSDSSDSSFSRSPPPGKQDSSDDVRRVQRREKNRIAAQKSRQRQTQKADTLHLESED
LEKQNAALRKEIKQLTEELKYFTSVLNSHEPLCSVLAASTPSPPEVVYSAHAFHQPHVSSPRFQP

纯度:> 95% by HPLC

浓度:1 mg/ml (determined by BCA assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.