

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

BCAR1, 465-848aa, Human, His tag, E.coli

产品货号: TP01301

第三版

别名:breast cancer anti-estrogen resistance 1, CAS, CAS1, CASS1, CRKAS, P130Cas

描述:BCAR1, or CAS, is an Src family kinase substrate involved in various cellular events, including migration, survival, transformation, and invasion. Recombinant human BCAR1 protein, fused to His-tag at N-terminus, was expressed in E.coli.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.4M UREA, 10% glycerol

分子量:43.9kDa (407aa)

序列:

MGSSHHHHHHSSGLVPRGSHMGSRLQQGVSAITVAHLLDLAGSAGATGSRWSPSEPQEPLVQDLQAAVA AVQSAVHELL
EFARSAVGNAAHTSDRALHAKLSRQLQKMEDVHQTLVAHGQALDAGRGGSGATLEDLDRLVACSRVAPEDAKQLASFL
HGNASLLFRRTKATAPGPEGGTLHPNPTDKTSSIQSRPLPSPPKFTSQDSPDGQYENSEGGWMEDYDYVHLQGKEEFKT
QKELLEKGSITRQGKSQLELQQLKQFERLEQEVSRPIDHDLANWTPAQPLAPGRTGGLGPSDRQLLLFYLEQCEANLTTLT
NAVDAFFTAVATNQPPKIFVAHSKFVILSAHKLFIGDTLRQAKAADVRSQVTHYSNLLCDLLRGIVATTKAAALQYSPSPS
AAQDMVE

纯度:> 95% by HPLC

浓度:0.5 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.