

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

BCL2L11, 1-138 aa Human, His tag, E.coli

产品货号: TP01313

第三版

别名: BCL2-like 11, BAM, BIM, BIM-alpha6, BIM-beta6, BIM-beta7, BimEL, BimL, BOD.

描述: BCL2L11 (Bcl-2-like protein 11) belongs to the Bcl-2 family. BCL-2 family members form hetero- or homodimers and act as anti- or pro-apoptotic regulators that are involved in a wide variety of cellular activities. The expression of BCL2L11 can be induced by nerve growth factor (NGF), as well as by the forkhead transcription factor FKHR-L1, which suggests a role of this gene in neuronal and lymphocyte apoptosis.

配方: Liquid. In 20 mM Tris-HCl buffer, pH8.0, 2M Urea, 20% glycerol 5mM DTT, 300mM NaCl

分子量: 18.5kDa (162aa)

序列:

MGSSHHHHHHSSGLVPRGSHMGSHMAKQPSDVSSECDREGRQLQPAERPPQLRPGAPTSLQTEPQDRSPAPMSCDKSTQ
TPSPPCQAFNHYLSAMASMRQAEPADMRPEIWIAQELRRIGDEFNAYYARRVFLNNYQAAEDHPRMVILRLLRYIVRLVW
RMH

纯度: > 95% by HPLC

浓度: 0.5mg/ml

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.