

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

BCL2L2, 1-172aa, Human, His tag, E.coli

产品货号: TP01314

第三版

别名:BCL2-like protein 2, Apoptosis regulator Bcl-W, BCLW.

描述:BCL2-like protein 2, also known as BCL2L2, is a member of the Bcl-2 family of proteins that regulates outer mitochondrial membrane permeability. BCL2L2 is an anti apoptotic member that prevents release of cytochrome c from the mitochondria intermembrane space into the cytosol.

配方:Liquid. In 20mM Tris-HCl buffer (pH7.5) containing 100mM NaCl, 10% glycerol

分子量:20.9 kDa (192aa)

序列:

MGSSHHHHHHSSGLVPRGSHMATPASAPDTRALVADFVGYKLRQKGYVCGAGPGEGPAADPLHQAMRAAGDEFETRF
RRTFSDLAAQLHVTPGSAQQRFTQVSEDLFQGGPNWGRLVAFFVFGAALCAESV NKEMEPLVGQVQEW MVAYLETRLA
DWIHSSGGWAEFTALYGDGALEEARRLREGNWASVRT

纯度:> 95% by HPLC

浓度:1 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储:+4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.