

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

BLyS, 134-285 aa Human, His-tagged, Recombinant, E.coli

产品货号: TP01346

第三版

别名: B lymphocyte stimulator, B-cell activating factor (BAFF), CD257, TALL1, THANK, ZTNF4, TNFSF13B, TNFSF20.

描述: B lymphocyte stimulator (BLyS), a member of the tumor necrosis factor (TNF) superfamily of cytokines, induces B-cell proliferation and immunoglobulin secretion and is a key regulator of peripheral B-cell populations in vivo. BLyS is made by immune-cells called monocytes and macrophages.

配方: Liquid. In 20 mM Tris-HCl buffer (pH 8.0) containing 5 mM DTT

分子量: 21 kDa (190 aa)

序列:

MRGSHHHHHHGMASMTGGQQMGRDLYDDDDKDRWGSHMAVQGPEETVTQDCLQLIADSETPTIQKGSYTFVPWLL
SFKRGSAL EEKENKILVKETGYFFIYGQVLYTDKTYAMGHLIQRKKVHVFGDEL SLVTLFRCIQNMPETLPNNSCYSAGIAK
LEEGDELQLAIPRENAQISLDGDVTFFGALKLL

纯度: > 95% by HPLC

浓度: 1 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.