

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

BMF, 1-129aa, Human, T7-tagged, Recombinant, E.coli

产品货号: TP01347

第三版

别名: Bcl2 modifying factor, isoform 3

描述: Bcl2 modifying factor, also known as BMF belongs to the Bcl2 protein family of apoptosis mediators. Bmf is constitutively expressed in many tissues. This protein contains a single Bcl2 homology domain 3 (BH3), and has been shown to bind Bcl2 proteins and function as an apoptotic activator.

配方: Liquid. In 20 mM Tris-HCl buffer (pH 7.5) containing, 1mM DTT, 10% glycerol

分子量: 15.6 kDa (144aa)

序列:

MASMTGGQQMGRGSHMEPSQCVEELEDVDFQPEDGEPVTQPGSLLSADLFAQSLLDCPLSRLQLFPLTHCCGPGLRPTSQ  
EDKATQTLSPASPSQGVMLPCGVTEEPQRLFYAPAEPKSCVVADPPLPAQPCFEWRREQERGRP

纯度: > 95% by HPLC

浓度: 1 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.