



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

BMP5, 317-454 aa, Human, Recombinant, E.coli

产品货号: TP01351

第三版

别名: Bone morphogenetic protein 5

描述: BMP5 is a member of the bone morphogenetic protein family which is part of the transforming growth factor-beta superfamily. The superfamily includes large families of growth and differentiation factors. BMP5 may play a role in certain cancers.

配方: Liquid. In 10mM Sodium Citrate buffer (pH3.5) containing 10% Glycerol.

分子量: 15.7 kDa (139aa)

序列:

MAANKRKNQNRNKSSSHQDSSRMSSVGDYNTSEQKQACKKHELYVSFRDLGWQDWIIAPEGYAAFYCDGEC SFPLNAH  
MNATNHAI VQTLVHLMFPDHVPKPCCAPTKLNAISVLYFDDSSNVILKKYRNMVVRSCGCH

纯度: > 95% by HPLC

浓度: 1 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.