

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

BOLA3, 1-107aa, Human, His tag, E.coli

产品货号: TP01359

第三版

别名: BolA-like protein 3 isoform1 , MMDS2

描述: BOLA3 also known as BolA-like protein 3 isoform1. BOLA3 is novel non-classical secreted protein and the secretion was not dependent on its predicted signal peptide. Human BolAs belongs to three different groups with functional divergence of BolA subfamily, where the different helix-turn-helix motif among BOLA 1, BOLA 2, and BOLA3 could be responsible for their functional divergence. Recombinant human BOLA3 was expressed in E.coli and purified by using conventional chromatography techniques.

配方: Liquid. In Phosphate buffered saline (pH7.4) containing 10% glycerol

分子量: 14.5kDa (130aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSMMAAWSPAAAAPLLRGIRGLPLHHRMFATQTEGELRVVTQILKEKFP RATAIKVTDISG  
GCGAMY EIKIESEEFKEKRTVQQHQMVNQALKEEIKEMHGLRIFTSVPKR

纯度: > 95% by HPLC

浓度: 1mg/ml (determined by Absorbance at 280nm)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.