

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

BRD1, 556-688aa Human, His tag, E.coli

产品货号: TP01365

第三版

别名: Bromodomain-containing protein 1, BRL, BRPF1, BRPF2.

描述: BRD contains a bromodomain, a sequence motif often found in transcriptional coactivators, and localizes to the nucleus in testis and several other cell types. Recombinant human BRD1 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方: Liquid. In Phosphate buffer saline (pH 7.4) containing 10% glycerol 1mM DTT

分子量: 17.8 kDa (156aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSEQVAMELRLTPLTVLLRSVLDQLQDKDPAIFAQPVSLKEVPDYLDHIKHPMDFAT
MRKRLEAQGYKNLHEFEEDFDLIIDNCMKYNARDTVFYRAAVRLRDQGGVVLQRARREVD SIGLEEASGMHLPERPA

纯度: > 95% by HPLC

浓度: 0.5 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.