



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

BTF3L4, 1-78aa Human, His tag, E.coli

产品货号: TP01375

第三版

别名: Transcription factor BTF3 homolog 4 isoform 3, Basic transcription factor 3-like 4

描述: BTF3L4 belongs to the NAC-beta family and contains 1 NAC-A/B (NAC-alpha/beta) domain. Recombinant human BTF3L4 protein isoform 3, fused to His-tag at N-terminus, was expressed in E.coli .

配方: Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 0.4M Urea

分子量: 11.3 kDa (101aa)

序列:

MGSSHHHHHSSGLVPRGSHMGSMNQEKLAKLQAQVRIGGKGTARRKKKVVHRTATADDKQLQSSLKKLAVNNIAGI
EESWTVKHQNQKTLMRKMMMFQIL

纯度: > 95% by HPLC

浓度: 1.0 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.