



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

CCDC25, 1-208 aa, Human, His tag, E.coli

产品货号: TP01446

第三版

别名:CCDC25, coiled-coil domain containing 25

描述:CCDC25 is a 208 amino acid protein that belongs to the CCDC25 family. Recombinant human CCDC25 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% glycerol, 1mM DTT

分子量:26.9 kDa (231aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSMVFYFTSSSVNSSAYTIYMGKDKYENEDLIKHWGPEDIWFHVDKLSSAHVYLRRLHK
GENIEDIPKEVLMDCAHLVKANSIQGCKMNNVNVVYTPWSNLKKTADMVVGQIGFHRQKDVKIVTVEKKVNEILNRLE
KTKVERFPDLAAEKECRDREERNEKKAQIQEMKKREKEEMKKKREMDLRSYSSLMKVENMSSNQDGNDSDEFM

纯度:> 95% by HPLC

浓度:1 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.