

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

CCL1, 24-96aa, Human, His tag, E.coli

产品货号: TP01450

第三版

别名: Chemokine C-C motif ligand 1

描述: CCL1 (chemokine C-C motif ligand 1) is a small glycoprotein secreted by activated T cells that belongs to a family inflammatory cytokines known as chemokines. CCL1 attracts monocytes, NK cells, and immature B cells and dendritic cells by interacting with a cell surface chemokine receptor called CCR8.

配方: Liquid. In 20mM Tris-HCl buffer (pH 7.5) containing 2mM DTT, 10% glycerol, 50mM NaCl

分子量: 10.8 kDa (94aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMKSMQVPFSRCCFSFAEQEIPLRILCYRNTSSICSNGLIFKLKRGKEACALDTVGVVQR  
HRKMLRHCPKSRK

纯度: > 95% by HPLC

浓度: 0.5 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.