



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

CCL13, 24-98aa, Human, His tag, E.coli

产品货号: TP01451

第三版

别名:C-C motif chemokine 13, CKb10, MCP-4, NCC-1, NCC1, SCYA1, SCYL1, CK-beta-10, Small-inducible cytokine A13

描述:CCL13 is a small cytokine belonging to the CC chemokine family that is also known as thymus and activation regulated chemokine (TARC). This protein induces chemotaxis in monocytes, eosinophils, T lymphocytes, and basophils by binding cell surface G-protein linked chemokine receptors such as CCR2, CCR3 and CCR5.

配方:Liquid. In Phosphate Buffered Saline pH7.4 containing 10% glycerol, 0.1M NaCl,1mM DTT

分子量:10.8 kDa (96aa), confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMQPDALNVPSTCCFTFSSKKISLQRLKSYVITTSRCPQKAVIFRTKLGKEICADPKEKWWVQ  
NYMKHLGRKAHTLKT

纯度:> 95% by HPLC

浓度:0.25 mg/ml (determined by BCA assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.