

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

CCL15, 22-113aa, Human, E.coli

产品货号: TP01453

第三版

别名: Chemokine (C-C motif) ligand 15, HCC-2, HMRP-2B, Lkn-1, LKN1, MIP-1d, MIP-5NCC-3, NCC3, SCYA15, SCYL3, SY15

描述: Chemokine (C-C motif) ligand 15, also known as CCL15, is a small cytokine belonging to the CC chemokine family that is also known as thymus and activation regulated chemokine. CCL15 induces changes in intracellular calcium concentration in monocytes and is thought to act through the CCR1 receptor.

配方: Liquid. In Phosphate buffered saline (pH7.4), 10% glycerol

分子量: 10.3kDa (93aa) confirmed by MALDI-TOF

序列:

MQFINDAETELMMSKLPLENPVVLNSFHFAADCCTSYISQSIPCSLMKSYFETSSECSKPGVIFLTKKGRQVCAKPSGPGVQ  
DCMKKLPYSI

纯度: > 95% by HPLC

浓度: 0.5 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.