

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

CCL17, 24-94 aa, Human, His-tagged, Recombinant, E.coli

产品货号: TP01454

第三版

别名: Chemokine (C-C motif) ligand 17, A-152E5.3, ABCD-2, SCYA17, TARC

描述: CCL17 is a small cytokine belonging to the CC chemokine family that is also known as thymus and activation regulated chemokine (TARC). This protein is involved in immunoregulatory and inflammatory processes. It is expressed constitutively in thymus, but only transiently in phytohemagglutinin-stimulated peripheral blood mononuclear cells.

配方: Liquid. In Phosphate Buffered Saline (pH7.4) containing 10% glycerol

分子量: 10.3 kDa (92aa)

序列:

MGSSHHHHHHSSGLVPRGSHMARGTNVGRECCLEYFKGAIPLRKLKTWYQTSEDCSRDAIVFVTVQGGAICSDPNNKR
VKNAV KYLQSLERS

纯度: > 95% by HPLC

浓度: > 95% by SDS - PAGE

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.