



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

CCL18, 22-89aa, Human, His tag, E.coli

产品货号: TP01455

第三版

别名:C-C motif chemokine 18, AMAC-1, AMAC1, CKb7, DC-CK1, DCCK1, MIP-4, PARC, SCYA18

描述:CCL18, also known as chemokine (C-C motif) ligand 18, a chemotactic factor that attracts lymphocytes but not monocytes or granulocytes may be involved in B cell migration into B cell follicles in lymph nodes. It attracts naive T lymphocytes toward dendritic cells and activated macrophages in lymph nodes, has chemotactic activity for naive T cells, CD4+ and CD8+ T cells and thus may play a role in both humoral and cell-mediated immunity responses.

配方:Liquid. In 10mM Sodium Citrate pH3.5, 10% Glycerol

分子量:10.4 kDa (93aa), confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSHMQVGTNKELCCLVYTSWQIPQKFIVDYSETSPQC PKPGVILLTKRGRQICADPNK
KWVQKYISDLKLNA

纯度:> 95% by HPLC

浓度:0.25 mg/ml (determined by BCA assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.