

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

CCL19, 22-98aa, Human, T7 tag, E.coli

产品货号: TP01456

第三版

别名:C-C motif chemokine 19, CKb11, ELC, MIP-3b, MIP3B, SCYA19

描述:CCL19 is a small cytokine belonging to the CC chemokine family that is also known as thymus and activation regulated chemokine (TARC). This protein is involved in immunoregulatory and inflammatory processes. It elicits its effects on its target cells by binding to the chemokine receptor chemokine receptor CCR7.

配方:Liquid. In Phosphate Buffered Saline pH7.4 containing 10% glycerol

分子量:10.4 kDa (93aa) confirmed by MALDI-TOF

序列:

MASMTGGQQMGRGSHMGTNDAEDCCLSVTQKPIPGYIVRNHFHYLLIKDGCRVPAVVFTTLRGRQLCAPPDQPWVERII  
QRLQRRTSAKMKRRSS

纯度:> 95% by HPLC

浓度:0.5 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.