

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

CCL20, 27-96aa, Human, His tag, E.coli

产品货号: TP01458

第三版

别名:C-C motif chemokine 20, CKb4, LARC, MIP-3a, MIP3A, SCYA20, ST38

描述:CCL20, also known as Macrophage Inflammatory Protein-3 (MIP3A), is a small cytokine belonging to the CC chemokine family. CCL20 protein is strongly chemotactic for lymphocytes and weakly attracts neutrophils. It is implicated in the formation and function of mucosal lymphoid tissues via chemoattraction of lymphocytes and dendritic cells towards the epithelial cells surrounding these tissues.

配方:Liquid. In Phosphate Buffered Saline pH7.4 containing 10% glycerol

分子量:10.3 kDa (91aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMASNFDCCCLGYTDRILHPKFIVGFTRQLANEGCDINAIIFHTKKKLSVCANPKQTWVKYI
VRLLSKKVKNM

纯度:> 95% by HPLC

浓度:0.25 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.