



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

CCL21, 24-134aa, Human, His tag, E.coli

产品货号: TP01459

第三版

别名:C-C motif chemokine 21, 6Ckine, CKb9, ECL, SCYA21, SLC, TCA4.

描述:CCL21 is member of the chemokine superfamily and the subfamily of CC chemokines that has an aspartate-cysteine-cysteine-leucine motif near its amino-terminus. This protein is involved in inhibits hemopoiesis and stimulates chemotaxis. CCL21 is chemotactic in vitro for thymocytes and activated T-cells, but not for B-cells, macrophages, or neutrophils.

配方:Liquid. In 20mM Tris buffer(pH 8.0) containing 10% glycerol 1mM DTT.I

分子量:14.5 kDa (132aa)

序列:

MGSSHHHHHHSSGLVPRGSHMSDGGQAQDCCLKYSQRKIPAKVVRYSYRKQEPSLGC SIPAILFLPRKRSQAELCADPKELW  
VQQLMQHLDKTPSPQKPAQGCRKDRGASKTGKKKGKSGKCKRTERSQTPKGP

纯度:> 95% by HPLC

浓度:0.5 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.