



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

CCL22/MDC, 25-93 aa, Human, His-tagged, Recombinant, E.coli

产品货号: TP01460

第三版

别名: Chemokine (C-C motif) ligand 22 / Macrophage-Derived Chemokine, SCYA22, STCP-1

描述: CCL22/MDC is a small cytokine belonging to the CC chemokine family. This protein displays chemotactic activity for natural killer cells, chronically activated T lymphocytes, monocytes and dendritic cells. It interacts with cell surface chemokine receptors CCR4.

配方: Liquid. In phosphate-buffered Saline(PBS) pH7.4, 10% glycerol

分子量: 10.3 kDa (90aa)

序列:

MGSSHHHHHSSGLVPRGSHMGYPYGANMEDSVCCRDYVRYRLPLRVVKHFYWTSDSCPRPGVLLTFRDKEICADPRV  
PWVKMILNKLSQ

纯度: > 95% by HPLC

浓度: 0.5 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.