

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

CCL27, 25-112aa, Human, E.coli

产品货号: TP01463

第三版

别名:C-C motif chemokine 27, ALP, CTACK, CTAK, ESKINE, ILC, PESKY, SCYA27

描述:Cutaneous T cell-attracting chemokine CCL27 is a member of the chemokine superfamily and the subfamily of r C-C chemokines that binds chemokine receptor CCR10. Chemokines are a superfamily of small secreted proteins that attract their targets by interacting with G protein-coupled receptors expressed on the migrating cell.

配方:Liquid. In 10 mM Sodium Citrate (pH 3.5), containing 10% Glycerol.

分子量:10.3 kDa (89aa), confirmed by MALDI-TOF

序列:

MFLLPPSTACCTQLYRKPLSDKLLRKVIQVELQEADGDCHLQAFVLHLAQRSICIHQPQNPSLSQWFEHQERKLHGTLPKLN
FGMLRKMG

纯度:> 95% by HPLC

浓度:0.5 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.