



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

CCL28, 23-127 aa, Human, His-tagged, Recombinant, E.coli

产品货号: TP01464

第三版

别名: Chemokine (C-C motif) ligand 28, CCK1, MEC, SCYA28

描述: CCL28, also known as chemokine (C-C motif) ligand 28 or mucosae-associated epithelial chemokine (MEC), belongs to the subfamily of small cytokine CC proteins. CCL28 is expressed by columnar epithelial cells in the gut, lung, breast and the salivary glands and drives the mucosal homing of T and B lymphocytes that express CCR10, and the migration of eosinophils expressing CCR3.

配方: Liquid. In 10mM Sodium Citrate pH3.5 containing 10% Glycerol

分子量: 14.3 kDa (126aa)

序列:

MGSSHHHHHHSSGLVPRGSHMILPIASSCCTEVSHHISRLLERVNMCRIQRADGDCDLAAVILHVKRRRICVSPHNHTV  
KQWMKVQAAKNGKGNVCHRKKHHGKRNSNRAHQGKHETYGHKTPY

纯度: > 95% by HPLC

浓度: 1 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.