

本公司提供的电版本说明书仅供参考，实验请以收到的纸质手册为准。

CCL3, 24-92aa, Human, His tag, E.coli

产品货号: TP01465

第三版

别名:C-C motif chemokine 3, G0S19-1, LD78ALPHA. MIP-1-alpha. MIP1A. SCYA3.

描述:CCL3, also known as, macrophage inflammatory protein-1alpha (MIP-1alpha), is a cytokine belonging to the CC chemokine family that is involved in the acute inflammatory state in the recruitment and activation of polymorphonuclear leukocytes.

配方:Liquid. In 10mM Sodium Citrate (pH3.5) containing 10% glycerol

分子量:10.0 kDa (90aa), confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMSLAADTPTACCFSYTSRQIPQNFADYFETSSQCCKPGVIFLTKRSRQVCADPSEEWVQK
YVSDLELSA

纯度:> 95% by HPLC

浓度:0.25 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.