

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

CCL3, 24-92aa, Mouse, His tag, E.coli

产品货号: TP01466

第三版

别名:C-C motif chemokine 3, AI323804, G0S19-1, LD78alpha, MIP-1alpha, MIP1-(a), MIP1-alpha, Mip1a, Scya3.

描述:C-C motif chemokine 3, also known as CCL3, belongs to the intercrine beta (chemokine CC) family. CCL3 has a potent chemotactic activity for eosinophils. This protein is binding to a high-affinity receptor activates calcium release in neutrophils.

配方:Liquid. In 20mM Sodium citrate buffer (pH 3.5) containing 10% glycerol

分子量:10.4 kDa (93aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHSSGLVPRGSHMGSMAPYGADTPTACCFYSRKIPRQFIVDYFETSSLCSQPGVIFLTKRNRQICADSKETW  
VQEYITDLELNA

纯度:> 95% by HPLC

浓度:0.5 mg/ml

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.