

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

CCL3L1, 26-93aa, Human, His tag, E.coli

产品货号: TP01467

第三版

别名:C-C motif chemokine 3-like 1, CCL3L3, D17S1718, G0S19-2, LD78, LD78BETA, MIP1AP, SCYA3L; SCYA3L1

描述:CCL3L1 is a small cytokine belonging to the CC chemokine. This protein is involved in immunoregulatory and inflammatory processes. It binds to several chemokine receptors including chemokine binding protein 2 (CCBP2 or D6) and chemokine (C-C motif) receptor 5 (CCR5).

配方:Liquid. In Phosphate Buffered Saline pH7.4 containing 10% glycerol

分子量:10.0 kDa (90aa), confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMSLAADTPTACCFSYTSRQIPQNFADYFETSSQCCKPSVIFLTKRGRQVCADPSEEWVQK
YVSDLELSA

纯度:> 95% by HPLC

浓度:0.25 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.