

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

CCL4L1, 24-92aa, Human, His tag, E.coli

产品货号: TP01468

第三版

别名:C-C motif chemokine 4-like, AT744.2, CCL4L, CCL4L2, LAG-1, LAG1, SCYA4L

描述:CCL4L1 (C-C motif chemokine 4-like) belongs to intercrine beta (chemokine CC) family. This protein is similar to CCL4 which inhibits HIV replication in peripheral blood monocytes that express CCR5. Recombinant human CCL4L1 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by conventional chromatography, after refolding of the isolated inclusion bodies in a renaturation buffer.

配方:Liquid. In 10mM sodium citrate (pH 3.5), 10% glycerol

分子量:10.5kDa (94aa) confirmed by MALDI-TOF (Molecular weight on SDS-PAGE will appear higher)

序列:

MGSSHHHHHHSSGLVPRGSHMGSHMAPMGSDPPTACCFSTARKLPRNFVVDYETSSLCSQPAVVFQTKRGKQVCAD
PSESWVQEYVYDLELN

纯度:> 95% by HPLC

浓度:0.5mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.