

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

CD1E, 32- 305aa Human, His tag, E.coli

产品货号: TP01483

第三版

别名:T-cell surface glycoprotein CD1e, , CD1A, R2

描述:CD1E is a member of the CD1 family of transmembrane glycoproteins, which are structurally related to the major histocompatibility complex (MHC) proteins and form heterodimers with beta-2-microglobulin. The CD1 proteins mediate the presentation of primarily lipid and glycolipid antigens of self or microbial origin to T cells. The human genome contains five CD1 family genes organized in a cluster on chromosome 1. The CD1 family members are thought to differ in their cellular localization and specificity for particular lipid ligands. CD1E localizes within Golgi compartments, endosomes, and lysosomes, and is cleaved into a stable soluble form. The soluble form is required for the intracellular processing of some glycolipids into a form that can be presented by other CD1 family members. Recombinant human CD1E protein, fused to His-tag at N-terminus, was expressed in E.coli.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 0.4M Urea

分子量:33.1 kDa (297aa)

序列:

MGSSHHHHHSSGLVPRGSHMGSEEQLSFRMLQTSSFANHSSWAHSESGWLGDLQTHGWDTVLGTIRFLKPWSHG NFS
KQELKNLQSLFQLYFHFSFIQIVQASAGQFQLEYPFEIQILAGCRMNAPQIFLN MAYQGSDFLSFQGISWEPSPGAGIRAQNIC
KVLNRYLDIKEILQSLLGHTCPRFLAGLMEAGESELKRKVKPEAWLSCGPSPGPGRLQLVCHVSGFYKPKVWVMWMRGE
QEQRGTQRGDVLPNADETWYLRATLDVAAGEAAGLSCRVKHSSLGGHDLIIHWGGYS

纯度:> 95% by HPLC

浓度:1mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.