

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

CD2, 25-209aa, Human, His tag, E.coli

产品货号: TP01484

第三版

别名:T-cell surface antigen CD2, SRBC, T11

描述:CD2, also known as T-cell surface antigen T11/Leu-5, is a cell adhesion molecule found on the surface of T cells and natural killer (NK) cells. It interacts with other adhesion molecules, such as lymphocyte function-associated antigen-3 (LFA-3/CD58) in humans, or CD48 in rodents, which are expressed on the surfaces of other cells. In addition to its adhesive properties, CD2 also acts as a co-stimulatory molecule on T and NK cells. CD2 is a specific marker for T cells and NK cells, and can therefore be used in immunohistochemistry to identify the presence of such cells in tissue sections. The great majority of T cell lymphomas and leukaemias also express CD2, making it possible to use the presence of the antigen to distinguish these conditions from B cell neoplasms. Recombinant human CD2 protein, fused to His-tag at N-terminus, was expressed in E.coli.

配方:Liquid. 20mM Tris-HCl buffer (pH8.0) containing 10% glycerol 0.4M Urea

分子量:23.8 kDa (209aa)

序列:

MGSSHHHHHSSGLVPRGSHMGSMKEITNALETWGALGQDINLDIPSFQMSDDIDDIKWEKTSKDKKIAQFRKEKETFK
EKDITYKLFKNGTLKIKHLKTDDQDIYKVSIDYDTKGKNVLEKIFDLKIQERVSKPKISWTCINTTLTCEVMNGTDPELNLYQ
DGKHLKLSQRVITHKWTTSLSAFKCTAGNKVSKESSVEPVSCPEKGLD

纯度:> 95% by HPLC

浓度:1 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.