

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

CD244, 19-224 aa Human, His tag, E.coli

产品货号: TP01493

第三版

别名:CD244 molecule, natural killer cell receptor 2B4, 2B4, NAIL, NKR2B4, Nmrk, SLAMF4.

描述:CD244 (Natural Killer Cell Receptor 2B4), also known as Cluster of Differentiation 244, contains 2 Iglike (immunoglobulin-like) domains. A role for the subtypes of CD2 Ig superfamily receptors has been recently demonstrated in eosinophilic inflammation in experimental asthma and atopic asthmatics.

配方:Liquid. In 20 mM Tris-HCl buffer, pH8.0, 0.4M Urea

分子量:25.5kDa (230aa)

序列:

MGSSHHHHHSSGLVPRGSHMGSHGKGCQGSADHVVSISGVPLQLQPNSIQTKVDSIAWKLLPSQNGFHHILKWENG  
SLPSNTSNDRFSEFIVKNLSLLIKAAQQQDSGLYCLEVTSISGKVQTATFQVFVFDKVEKPRLQGQGGKILDRGRCQVALSCLV  
SRDGNVSYAWYRGSKLIQTAGNLTYLDEEVDINGTHTYTCNVSNPVSWEHTLNLTQDCQNAHQEFRFWP

纯度:> 95% by HPLC

浓度:0.5mg/ml

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.