

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

CD27, 21-191aa, Human, His tag, E.coli

产品货号: TP01497

第三版

别名:CD27 antigen, S152, T14, TNFRSF7, Tp55.

描述:TNFRSF7, also known as CD27, is a member of the TNF-receptor superfamily. This receptor is required for generation and long-term maintenance of T cell immunity. It binds to ligand CD70, and plays a key role in regulating B-cell activation and immunoglobulin synthesis.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% glycerol

分子量:21.8 kDa (196aa)

序列:

MGSSHHHHHHSSGLVPRGSHMGSHMTPAPKSCPERHYWAQGKLCCQMCEPGTFLVKDCDQHRKAAQCDPCIPGVSF
SPDHHTRPHCESCRHCNSGLLVRNCTITANAECACRNGWQCRDKECTECDPLPNPSLTARSSQALSPHPQPTHLPYVSE
MLEARTAGHMQTLADFRQLPARTLSTHWPPQRSRLCSDSFIR

纯度:> 95% by HPLC

浓度:1 mg/ml

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.