



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

CD274, 19-239aa, Mouse, His tag, E.coli

产品货号: TP01499

第三版

别名:CD274 antigen, A530045L16Rik, B7h1, Pcd11l1, Pcd1lg1,Pdl1.

描述:CD274, also known as programmed cell death 1 ligand 1 (PDCD1L1), is involved in the costimulatory signal, essential for T-cell proliferation and production of IL10 and IFNG, in an IL2-dependent and a PDCD1- independent manner.

配方:Liquid. 20mM Tris-HCl buffer (pH8.0) containing 10% glycerol, 0.15M NaCl

分子量:27.6kDa(246aa), confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSHMFTITAPKDLYVVEYGSNVTMECRFPVERELDLLALVVYWEKEDEQVIQFVAGE  
EDLKPQHSNFRGRASLPKDQLLKGNAAALQITDVKLQDAGVYCCIISYGGADYKRITLKVNPYRKINQRISVDPATSEHELI  
CQAEGYPEAEVIWTNSDHQPVSGKRSVTTSRTEGMLLNVTSSLRVNATANDVFYCTFWRSQPGQNHTAELIPELPATHP  
PQNRTH

纯度:> 95% by HPLC

浓度:1 mg/ml

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.