

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

CD68, 22-319aa, Human, His tag, E.coli

产品货号: TP01525

第三版

别名:Macrosialin , GP110, LAMP4, SCARD1.

描述:Macrosialin, also known as CD68, is a member of the lysosomal/endosomal-associated membrane glycoprotein (LAMP) family. It is found in cytoplasmic granules and in the cytoplasm of various non-hematopoietic tissues including liver and kidney tubules and glomeruli.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 1M Urea, 10% glycerol, 0.1M NaCl

分子量:34.1 kDa (323aa)

序列:

MGSSHPLHHHSSGLVPRGSHMGSHMNDCPHKKSATLLPSFTVTPTVTESTGTTSHRTTKSHKTTTHRTTTTGTTSHGPTT
ATHNPTTTTSHGNVTVHPTSNSTATSQGPSTATHSPATTSHGNATVHPTSNSTATSPGFTSSAHPEPPPPSPSPSPTSKETIGD
YTWTNGSQPCVHLQAQIQIRVMYTTQGGGEAWGISVLNPNKTKVQGSCEGAHPHLLLSFPYGHLISFGFMQDLQKVV
YLSYMAVEYNVSFPHAAQWTFSAQNASLRDLQAPLGQSFSCSNSSIIILSPAVHLDLLSLRLQAAQLPHTGVFGQSFSCPSDRS

纯度:> 95% by HPLC

浓度:1 mg/ml

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.