

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

CD70, 39-193aa, Human, His tag, E.coli

产品货号: TP01531

第三版

别名:CD70 antigen, CD70L, CD70LG, TNFSF7

描述:CD70 also known as CD70 antigen is cytokine that binds to CD27. This protein plays a role in T-cell activation and induces the proliferation of costimulated T-cells and enhances the generation of cytolytic T-cells. Recombinant human CD70, fused to His-tag at N-terminus, was expressed in E.coli

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol

分子量:19.5kDa (178aa)

序列:

MGSSHHHHHHSSGLVPRGSHMGSQRFAQAQQQLPLESLGWDVAELQLNHTGPQQDPRLYWQGGPALGRSFLHGPELD
KGQLRIHRDGIYMVHIQVTLAICSSTTASRHHPTTLAVGICSPASRSISLLRSLFHQGCTIASQRLTPLARGDTLCTNLGTLL
PSRNTDETFFGVQWVRP

纯度:> 95% by HPLC

浓度:1mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.