

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

CD84, 22-225aa, Human, His tag, E.coli

产品货号: TP01538

第三版

别名:SLAM family member 5, hCD84, LY9B, mCD84, SLAMF5.

描述:CD84 is a membrane glycoprotein that is a member of the signaling lymphocyte activation molecule (SLAM) family. This family forms a subset of the larger CD2 cell-surface receptor Ig superfamily. The protein is a homophilic adhesion molecule that is expressed in numerous immune cells types and is involved in regulating receptor-mediated signaling in those cells.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% glycerol

分子量:25.4 kDa (229aa)

序列:

MGSSHHHHHHSSGLVPRGSHMGSHMKDSEIFTVNGILGESVTFPVNIQEPRQVKIIAWTSKTSVAYVTPGDSETAPVVTV
THRNYERIHALGPNYNLVISDLRMDAGDYKADINTQADPYTTTKRYNLQIYRRLGKPKITQSLMASVNSTCNVTLTCS
VEKEEKVNTYNWSPLGEEGNVLQIFQTPEDQELTYTCTAQNPVSNNSDSISARQLCADIAMGFRTHHTG

纯度:> 95% by HPLC

浓度:1 mg/ml

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.