

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Cellubrevin, 1-77aa, Human, Recombinant, E.coli

产品货号: TP01576

第三版

别名: Vesicle-associated membrane protein 3, VAMP 3, CEB, Synaptobrevin-3

描述: Cellubrevin, also known as VAMP 3, is present in recycling endosomes and endosome-derived vesicles. This protein has been implicated in recycling of transferrin receptors to the plasma membrane, secretion of alpha-granules in platelets, recycling of T-cell receptors to the immunological synapses, and membrane trafficking during cell migration.

配方: Liquid in 20mM Tris-HCl pH 7.5, 10% glycerol

分子量: 8.7 kDa (77aa)

序列:

MSTGPTAATGSNRRLQQTQNQVDEVVDIMRVNVDKVLERDQKLSLDDRADALQAGASQFETSAAKLKRKYWWKNC  
K

纯度: > 95% by HPLC

浓度: 1 mg/ml

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.