



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

CGREF1, 20-301aa, Human, His tag, E.coli

产品货号: TP01598

第三版

别名: Cell growth regulator with EF-hand domain 1, CGR11

描述: MGSSHHHHHH SSGLVPRGSH MGSAPKDGVT RPDSEVQHQL LPNPFQPGQE QLGLLQSYLK GLGRTEVQLE
HLSREQVLLY L FALHDYDQS GQLDGLELLS MLTAALAPGA ANSPTTNPVI LIVDKVLETQ DLNGDGLMTP
AELINFPQVA LRHVEPGEPL APSPQEPQAV GRQSLLAKSP LRQETQEAPG PREEAKGQVE ARRESLDPVQ
EPGGQAEADG DVPGRGEAE GQAEAKGDAP GPRGEAGGQA EAEGDAPGPR GEAGGQAEAR ENGEEAKELP
GETLESKNTQ NDFEVHIVQV ENDEI

配方: Liquid. Phosphate Buffered Saline (pH7.4)

分子量: 32.3 kDa(305aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSAPKDGVT RPDSEVQHQLLPNPFQPGQEQLGLLQSYLKGLGRTEVQLEHLSREQVLL
YLFALHDYDQSGQLDGLELLSMLTAALAPGAANSPTTNPVILIVDKVLETQDLNGDGLMTPAELINFPQVALRHVEPGEPL
APSPQEPQAVGRQSLLAKSPLRQETQEAPGPREEAKGQVEARRESLDPVQEPGGQAEADGDVPGRGEAEAGQAEAKGDA
PGPRGEAGGQAEAEAGDAPGPRGEAGGQAEARENGEEAKELPGETLESKNTQ NDFEVHIVQV ENDEI

纯度: > 95% by HPLC

浓度: 0.5 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.