

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

CHAC2, 1-184aa, Human, His tag, E.coli

产品货号: TP01599

---

第三版

**别名:** Cation transport regulator-like protein 2, ChaC, cation transport regulator homolog 2

**描述:** CHAC2, also known as cation transport regulator-like protein 2, belongs to the chaC family. Recombinant human CHAC2 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

**配方:** Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 10% glycerol, 1mM DTT

**分子量:** 23.3 kDa (207aa) confirmed by MALDI-TOF

**序列:**

**纯度:** > 95% by HPLC

**浓度:** 0.5 mg/ml (determined by Bradford assay)

**内毒素:** <1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.