

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

CHMP1B, 1-199aa, Human, His tag, E.coli

产品货号: TP01606

第三版

别名:Charged multivesicular body protein 1b, C10orf2, C18-ORF2, C18orf2, CHMP1.5, hVps46-2, Vps46-2, Vps46B.

描述:CHMP1B, also known as charged multivesicular body protein 1b, is probable peripherally associated component of the endosomal sorting required for transport complex III (ESCRT-III) which is involved in multivesicular bodies (MVBs) formation and sorting of endosomal cargo proteins into MVBs.

配方:Liquid. 20mM Tris-HCl buffer (pH8.0) containing 20% glycerol, 0.1M NaCl, 1mM DTT

分子量:24.7 kDa (223aa), confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSHMSNMEKHLFNLKFAAKELSRSAKKCDKEEKAEKAKIKKAIQKGNMEVARIHAEN  
AIRQKNQAVNFLRMSARVDAVAARVQTAVTMGKVTKSMAGVVKSM DATLKT MNLEKISALMDKFEHQFETLDVQTQ  
QMEDTMSSTTTLTTPQNQVDMLLQEMADEAGLDLNMELPQGQTGSVGT SVASAEQDELSQRLARLRDQV

纯度:> 95% by HPLC

浓度:0.25mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.