

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

CKS2, 1-79aa, Human, T7-tagged, Recombinant, E.coli

产品货号: TP01634

---

第三版

别名: CDC28 Protein kinase 2, CKSHS2

描述: CKS2, also known as CDC28 Protein kinase 2 associates with the catalytic subunit of cyclindependent kinases and has therefore been assumed to play a direct role in cell cycle regulation. The function of CKS2 in somatic mammalian cells is not well understood although it is required for the first metaphase/anaphase transition during the meiosis.

配方: Liquid in 20mM Tris pH7.5, 20% Glycerol

分子量: 11 kDa (94 aa)

序列:

纯度: > 95% by HPLC

浓度: 1 mg/ml

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.