

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

CLEC4E, 41-219aa, Human, His tag, E.coli

产品货号: TP01642

第三版

别名:C-type lectin domain family 4, member E, MINCLE, CLECSF9

描述:CLEC4E also known as C-type lectin domain family 4, member E. CLEC4E is a C-type lectin receptor expressed by a variety of antigen presenting cells including macrophages, neutrophils, dendritic cells and B cells; a variety of stimuli including stress are known to induce the expression of CLEC4E. Recombinant human CLEC4E, fused to His-tag at N-terminus, was expressed in E.coli.

配方:Liquid. In 20mM Tris 8.0 containing 10% glycerol

分子量:22.9kDa (200aa)

序列:

MGSSHHHHHHSSGLVPRGSHMRCVVTFRIFQTCDEKKFQLPENFTELSCYNYGSGSVKNCCPLNWEYFQSSCYFFSTDTIS
WALSLKNCSAMGAHLVVINSQEEQEFLSYKKPKMREFFIGLSDQVVEGQWQWVDGTPLTKLSFWDVGEPNNIATLEDC
ATMRDSSNPRQNWNDVTCFLNYFRICEMVGINPLNKGKSL

纯度:> 95% by HPLC

浓度:0.25mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.