

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

CNOT7, 1-248aa, Mouse, His tag, E.coli

产品货号: TP01664

第三版

别名:CCR4-Not transcription complex, subunit 7, CAF-1.

描述:CNOT7, also known as CCR4-Not transcription complex, subunit 7, is a ubiquitous transcription factor. It is a component of the CCR4 complex which functions as a general transcription regulation complex. CNOT7 also binds to an anti-proliferative protein, B-cell translocation protein 1 (BTG1), which negatively regulates cell proliferation.

配方:Liquid. 20mM Tris-HCl buffer (pH8.0) containing 20% glycerol, 1mM DTT

分子量:31.1 kDa (271aa), confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGMPAATVDHSQRICEVWACNLDEEMKKIRQVIRKYNVYVAMDTEFPGVVARPIGEF
RSNADYQYQLLR CNVDLLKIIQLGLTFMNEQGEYPPGTSTWQFNFKFNLTEDMYAQDSIELLTTSGIQFKKHEEEGIETQYF
AELLMTSGVVLCEGVKWLFSHSGYDFGYLIKILTNSNLPEEELDFEILRLFFPVIYDVKYLKMKCKNLKMFEDHIDDAKYC
GHLYGLGSGSSYVQNGTGNAYEEEEASKQS

纯度:> 95% by HPLC

浓度:1 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.