



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

CNTF(Ciliary neurotrophic factor), 1-200aa, Human, His-tagged, Recombinant, E.coli

产品货号: TP01670

第三版

别名:Ciliary neurotrophic factor

描述:CNTF (Ciliary neurotrophic factor) is a polypeptide hormone whose actions appear to be restricted to the nervous system where it promotes neurotransmitter synthesis and neurite outgrowth in certain neuronal populations. This protein is a potent survival factor for neurons and oligodendrocytes and may be relevant in reducing tissue destruction during inflammatory attacks.

配方:Liquid. In 20 mM Tris-HCl buffer (pH 8.5) containing 1 mM DTT, 30% Glycerol, 0.2M NaCl

分子量:25kDa (220aa), confirmed by MALDI-TOF.

序列:

MGSSHHHHHHSSGLVPRGSHMAFTEHSPLTPHRRDLCSRSIWLARKIRSDLTALTESYVKHQGLNKNINLDSADGMPVAS
TDQWSELTEAERLQENLQAYRTFHVLLARLLEDQQVHFTPTEGDFHQAIHTLLQVAAFAYQIEELMILLEYKIPRNEADG
MPINVDGGLFEKKLWGLKVLQELSQWTVRSIHDLRFISSHQTGIPARGSHYIANNKKM

纯度:> 95% by HPLC

浓度:1 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.