

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

CNTF, 1-200aa, Human, E.coli

产品货号: TP01671

第三版

别名:Ciliary neurotrophic factor.

描述:CNTF (Ciliary neurotrophic factor) is a polypeptide hormone whose actions appear to be restricted to the nervous system where it promotes neurotransmitter synthesis and neurite outgrowth in certain neuronal populations. This protein is a potent survival factor for neurons and oligodendrocytes and may be relevant in reducing tissue destruction during inflammatory attacks.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 10% glycerol

分子量:22.9 kDa (200aa), confirmed by MALDI-TOF

序列:

MAFTEHSPLTPHRRDLCSRSIWLARKIRSDLTALTESYVKHQGLNKNINLDSADGMPVASTDQWSELTEAERLQENLQAY
RTFHVLLARLLEDQQVHFTPTGDFHQAIHTLLLQVAAFAYQIEELMILLEYKIPRNEADGMPINVGDGGLFEKKLWGLK
VLQELSQWTVRSIHDLRFISSHQTGIPARGSHYIANNKKM

纯度:> 95% by HPLC

浓度:1 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.