

本公司提供的电版本说明书仅供参考，实验请以收到的纸质手册为准。

CNTFR, 23-342 aa, Human, His tag, E.coli

产品货号: TP01674

第三版

别名:Ciliary neurotrophic factor receptor subunit alpha, CNTFR-alpha, CNTFR Alpha, CNTF receptor

描述:CNTFR also known as Ciliary neurotrophic factor receptor subunit alpha. This protein binds to CNTF. The alpha subunit provides the receptor specificity. Recombinant human CNTFR protein, was expressed in E.coli

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol

分子量:38.1kDa (341aa)

序列:

MGSSHHHHHHSSGLVPRGSHMQRHSPQEAPHVQYERLGSDVTLPCGTANWDAAVTWRVNGTDLAPDLLNGSQLVLH
GLELGHSGLYACFHRDSWHLRHQVLLHVGLPPREPVLSRSTYYPKGFYCSWHLPTPTYIPNTFNVTVLHGSKIMVCEKD
PALKNRCHIRYMHLFSTIKYKVSISVSNALGHNATAITFDEFTIVKPDPPENVVVARPVPSNPRRLEVTWQTPSTWPDPEFPL
KFFLRYPRLILDQWQHVELSDGTAHTITDAYAGKEYIIQVAAKDNEIGTWSWVAHAHATPWTEPRHLTTEAQAETT
TSTTSSLAPPPTTKICDPGELGS

纯度:> 95% by HPLC

浓度:1mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.