

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

COMMD6, 1-85aa, Human, His tag, E.coli

产品货号: TP01680

第三版

别名: COMM domain-containing protein 6 isoform b, Acrg

描述: COMMD6 contains 1 COMM domain. This protein down-regulates activation of NF-kappa-B and inhibits TNF-induced NFKB1 activation. Recombinant human COMMD6 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography.

配方: Liquid. 20mM Tris-HCl buffer (pH7.0) containing 20% glycerol 0.15M NaCl, 1mM DTT

分子量: 12.0 kDa(108aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSMEASSEPLDAKSDVTNQLVDFQWKLGMVSSDTCRSLKYPYVAVMLKVADHSG  
QVKTKCFEMTIPQFQNFYRQFKEIAAVIETV

纯度: > 95% by HPLC

浓度: 1 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.