

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

COX4NB, 1-210aa, Human, His tag, E.coli

产品货号: TP01693

第三版

别名:Neighbor of COX4, NOC4.

描述:Neighbor of COX4, also known as COX4NB, belongs to the UPF0172 (NOC4) family and is found on chromosome 16, adjacent to the gene that encodes COX4. This protein is expressed in liver, pancreas, heart, lung, kidney, brain, skeletal muscle, and placenta.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 10% glycerol, 0.1M NaCl

分子量:25.9 kDa (230aa), confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMPGVKLTQAYCKMVLHGAKYPHCAVNGLLVAEKQKPRKEHLPLGGPGAHTLTFVD
CIPLFHGTLALAPMLEVALTLIDSWCKDHSYVIAGYYQANERVKDASPNQVAEKVASRIAEGFSDTALIMVDNTKFTMDC
VAPTIHVYEHENRWRCRDPHHDYCEDWPEAQRISASLLDSRSYETLVDFDNHLDDIRNDWTNPEINKAVLHLC

纯度:> 95% by HPLC

浓度:1 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.