



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

CRABP1, 1-137 aa, Human, Recombinant, E.coli

产品货号: TP01702

第三版

别名: Cellular retinoic acid binding protein 1, CRABP, CRABP-I, CRABPI, RBP5

描述: CRABP1, also known as cellular retinoic acid-binding protein 1, is a member of specific carrier proteins for members of the vitamin A family. This protein is assumed to play an important role in retinoic acid-mediated differentiation and proliferation processes.

配方: Liquid. In 20mM Tris-HCl buffer (pH8.0) containing 10% glycerol

分子量: 15.5 kDa (137 aa)

序列:

MPNFAGTWKMRSSNFDELLKALGVNAMLKRVAVAAASKPHVEIRQDGDQFYIKTSTTVRTTEINFKVGEGFEEETVDG  
RKCRSLATWENENKIHCTQTLLEGDGPITYWTRELANDELILTFGADDVVCTRIYVRE

纯度: > 95% by HPLC

浓度: 1.0 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.